









KunShan Adhesive Package International Co., Ltd






Add: 1504, No.9, Lane 231, Lvdi Road, Huaqiao Town, Kunshan, Jiangsu Province, China.





Tel:86 021 69781176 Fax:86 021 69781175





www.adhesive-pack.com Email:adhesive.package@gmail.com




Item No.	PICTURE	DESCRIPTION	FUNCTION	G.W. (kgs)	CTN SIZE (mm)	PACKING
7#		Description: Red Mixing Tip Mixing type: Dynamic Mixer Outlet Diameter:6.88mm Big entry diameter:7.22mm Small entry diameter:5.32mm Total Length: 49.22mm Ratio: 5:1	Mixing 3M heavy-body and precision dental impression materials. Suitable Machine:Pentamix™ Lite Automatic Mixing Unit	16.6	450*350*490	50pcs / bag 50bags / ctn 2500pcs / ctn
12# New		Description: Grey Mixing Tip Outlet Diameter: 8.78mm Big entry diameter:8.55mm Small entry diameter:5.75mm Total Length: 52.3mm Mixing type: Dynamic Mixer Ratio:5:1	Suitable Machine: Zhermack Modulmix Function: Mixing Heavy/Maxi Heavy/Maxi Monophase/Maxi Putty dental impression materials.	17	530*420*370	50pcs / bag 38bags / ctn 1900pcs / ctn
12#		Description: Yellow Mixing Tip Outlet Diameter:5.95mm Big entry diameter:8.35mm Small entry diameter:5.9mm Total Length: 52.3mm Mixing type: Dynamic Mixer Total Length: 55.3mm Ratio:5:1	Suitable Machine: Zhermack, Hereaus Kuzler, Dentsply, Conlte, Huge dental machines etc. Function: Mixing Regular/Heavy-body dental impression materials.	14	460*420*390	50pcs / bag 30bags / ctn 1500pcs / ctn


Item No.	PICTURE	DESCRIPTION	FUNCTION	G.W. (kgs)	CTN SIZE (mm)	PACKING
15#		Description: Blue Mixing Tip Outlet Diameter: 5.86mm Entry diameter:6.51mm Total Length:56.6mm Mixing type: Dynamic Mixer Ratio:1:1	Suitable Machine: Dentsply dental machine(old-type)	10.9	440*375*410	50pcs / bag 25bags / ctn 1250pcs / ctn
16#		Description:Green Mixing Tip Mixing type: Dynamic Mixer Outlet Diameter: 6mm Big entry diameter:6.9mm Small entry diameter:3.16mm Total Length:49.8mm Ratio:5:1	Suitable Machine: Coltene Whaledent dental machine(old-type)	10.9	440*375*410	50pcs / bag 25bags / ctn 1250pcs / ctn
17#		Description:White Mixing Tip Mixing type: Dynamic Mixer Outlet Diameter: 6mm Big entry diameter:6.9mm Small entry diameter:3.16mm Total Length:49.8mm Ratio:5:1	Suitable Machine: Coltene Whaledent dental machine(old-type)	10.9	440*375*410	50pcs / bag 25bags / ctn 1250pcs / ctn

Item No.	PICTURE	DESCRIPTION	FUNCTION	G.W. (kgs)	CTN SIZE (mm)	PACKING
19# (New-developed Product)		Description: DMG Mixing Tip Mixing type: Dynamic Mixer Outlet Diameter: 6mm Big entry diameter:6.7mm Small entry diameter:3.08mm Total Length: 52.5mm Ratio:5:1	Suitable Machine: DMG MixStar-eMotion	14	460*420*390	50pcs / bag 30bags / ctn 1500pcs / ctn
N1		Item no.: N-1 (yellow) Description: Dental Intra-oral Tip	Used together with 2# and 5# mixing tips, which is used for delivering dual-cured Dual-cured core build-up composite resin, light-body dental impression materials etc.	6.3	370*350*210	100pcs / bag 500bags / ctn 50000pcs / ctn
N2		Item no.: N-2 (clear) Description:Dental Intra-oral Tip	Used together with 3# and 4# mixing tips, which is used for delivering regular-body dental impression materials etc.	6.3	370*350*210	100pcs / bag 500bags / ctn 50000pcs / ctn
N3		Item no.: N-3 (clear) Description: Dental Intra-oral Tips	Used together with 11# and 14# mixing tips, which is used for delivering dental cold-whitening cement, dual-cured composite resins, dental luting and cement etc.	6	370*350*210	100pcs / bag 500bags / ctn 50000pcs / ctn
N5		Item no.: N-5 (clear) Description: Dental Intra-oral Tip	Used together with 3# mixing tip, which is used for bite registration.	6.8	370*350*210	100pcs / bag 500bags / ctn 50000pcs / ctn

Item No.	PICTURE	DESCRIPTION	FUNCTION	G.W. (kgs)	CTN SIZE (mm)	PACKING
N7		Item no.: N-7 (clear) Description: 3M Dental Application Tip	Used together with 3M capsule to delivery encapsulated glass ionomer products.	6.3	370*350*210	10pcs / bag 500bags / ctn
N8		Item no.: N-8 (clear) Description: Dental Intra-oral Tip	Used together with 11# and 14# mixing tips, which is used for delivering dual-cured radiopaque resin and dual-cored resin adhesive luting, adhesive resin cement	3.8	320*270*220	100pcs / bag 500bags / ctn 50000pcs / ctn
14ml 1:1 (New-developed Product)		Item no.: 14ml 1:1 Dual Syringes Material: PBT, PP Minimum wall thickness: \geq 1mm Part A Diameter : Outer-19mm Inner-17mm Part B Diameter : Outer-19mm Inner-17mm Length: 80mm	Packing various silicone impression materials in dental & beauty.	26	500*400*400	1000pcs / ctn
25ml 1:1		Item no.: 24ml 1:1 Dual Syringes Material: PBT, PP Minimum wall thickness: \geq 1mm Part A Diameter : Outer-19mm Inner-17mm Part B Diameter : Outer-19mm Inner-17mm Length: 80mm	Packing various silicone impression materials in dental & beauty.	20	500*400*400	800pcs / ctn

Item No.	PICTURE	DESCRIPTION	FUNCTION	G.W. (kgs)	CTN SIZE (mm)	PACKING
50ml 1:1		Size:50ml 1:1 Material:PP/PE Length:100mm Diameter(part A):22mm Diameter(part B) :22mm	Packing various dental impression materials, bite registration materials, soft/hard relining materials, dual-cured composite resins, self-curing temporary crown&bridge resins etc.	11	600*320*460	20pcs / bag 20bags / ctn
50ml 10:1		Size:50ml 10:1 Material:PA/PBT Length: 105mm Diameter(part A):30mm Diameter(part B):11mm	Packing dental temporary crown & bridge materials	17	460*420*410	20pcs / bag 20bags / ctn
CG50-1-1 II New Type For 1:1 / 2:1		Cartridge: two-component Capacity: 50ml (1:1,2:1) Material:Reinforced Nylon Frame+Plastic Rod	Dispensing any-brand dental impression materials, bite registration materials, soft/hard relining materials, dual-cured composite resins, self-curing temporary crown&bridge resins etc.	14.8	610*410*390	1pc/bag in box 50pcs/ctn
CG50-10-1 New Type For 10:1 / 4:1		Cartridge:two-component Capacity: 50ml (4:1,10:1) Material:Reinforced Nylon Frame+ Plastic Rod	Dispensing any-brand dental temporary crown & bridge materials.	14.3	610*410*390	1pc/bag in box 50pcs/ctn

Item No.	PICTURE	DESCRIPTION	FUNCTION	G.W. (kgs)	CTN SIZE (mm)	PACKING
Applicator/Applier		Item no.: 3M/GC Applierx Description: Made of Stainless Steel	Designed to be used with all encapsulated glass ionomer products to facilitate easy dispensing of 3M or GC's Amalgam mixture directly to tooth preparator.	16	Inner Box: 180*40*90 (0.15kg) Carton Size: 430*350*460	1pc/bag in box 100pcs/ctn
3M Dental Aplicap Capsule Activator		Item no.: 3M Capsule Activator Description: Made of Stainless Steel	Used for opening/activating the pluger of capsule. Insert the capsule into activator with the application tip point downward. Slowly depress the activating lever down to the top, and hold for 2 to 4 econds. It's essential to strictly comply with this proceduere to prevent premature curing of the material inside capsule.	16	Inner Box: 180*80*90 (0.3kg) Carton Size: 430*350*460	1pc/bag in box 50pcs/ctn
3M Set Shade Guide		Item No.:3M Set Shade Guide Description: 1.Made of light-cured composite resin 2.Standard Dental 3M shade guide with 8 colors (A1,A2,A3,A3.5,B1...) 3. For making resin teeth with differnt colors 4. Easy to compare	Used for making the resin teeth with differnt colors, so that the dentist can choose the teeth with right color for patient.	14	Inner Box: 165*80*30 Carton Size: 500*400*310	1pcs / box 120pcs / ctn

Item No.	PICTURE	DESCRIPTION	FUNCTION	G.W. (kgs)	CTN SIZE (mm)	PACKING
Amalgam Mixer		<p>Item No.: NSJ-DA110 Dental Amalgamator</p> <p>Features:</p> <ul style="list-style-type: none"> ·Silicone rubber buttons for comfortable touch ·LED display ·Two speed selections ·Noiseless and stable operation ·Easy operation and sturdy capsule holder ·Sealed mixing area with a safety cover ·Ergonomic design <p>Technical data:</p> <p>Input Voltage: 220V 110V 50HZ/60HZ</p> <p>Input power: 30W</p> <p>Mixing frequency: 2960RPM~4200RPM</p> <p>Mixing time: 4 to 40 seconds</p> <p>Dimensions: 295L X 250L X 230H</p> <p>Weight : 2.2KG</p>	<p>Used for mixing all mixed capsules including amalgam capsules. The universal high-speed amalgamator is controlled by a microprocessor, which ensures a precise, controllable and consistent mixing. With a safety cover, DA-pro will deactivate automatically when you open the front cover to access the capsule.</p>	9.5	<p>Inner Box: 260*270*190m m</p> <p>Carton: 570*300*420m m</p>	<p>1pcs / box 4pcs / ctn</p>